

Product datasheet for **GA111975**

Human tenomodulin (TNMD) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNMD CRISPRa kit - CRISPR gene activation of human tenomodulin
Symbol:	TNMD
UniProt ID:	<u>Q9H2S6</u>
Synonyms:	BRICD4; CHM1L; TEM
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)