

Product datasheet for **GA111983**

Human TMBIM1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMBIM1 CRISPRa kit – CRISPR gene activation of human transmembrane BAX inhibitor motif containing 1
Symbol:	TMBIM1
Components:	GA111983G1, TMBIM1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111983G2, TMBIM1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111983G3, TMBIM1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q969X1</u>
Synonyms:	LFG3; MST100; MSTP100; PP1201; RECS1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

