

Product datasheet for **GA111985**

Human **BIGM103 (SLC39A8)** activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC39A8 CRISPRa kit - CRISPR gene activation of human solute carrier family 39 member 8
Symbol:	SLC39A8
UniProt ID:	<u>Q9C0K1</u>
Synonyms:	BIGM103; CDG2N; LZT-Hs6; PP3105; ZIP8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)