


Product datasheet for **GA111986**

Human DUS1L activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DUS1L CRISPRa kit - CRISPR gene activation of human dihydrouridine synthase 1 like
Symbol:	DUS1L
Components:	GA111986G1, DUS1L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA111986G2, DUS1L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA111986G3, DUS1L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6P1R4</u>
Synonyms:	DUS1; PP3111
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)