

Product datasheet for GA111991

Human TINAGL1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TINAGL1 CRISPRa kit - CRISPR gene activation of human tubulointerstitial nephritis antigen like 1
Symbol:	TINAGL1
UniProt ID:	Q9GZM7
Synonyms:	ARG1; LCN7; LIECG3; TINAGRP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)