

Product datasheet for **GA112024**

Human MS4A6A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MS4A6A CRISPRa kit - CRISPR gene activation of human membrane spanning 4-domains A6A
Symbol:	MS4A6A
UniProt ID:	<u>Q9H2W1</u>
Synonyms:	4SPAN3; 4SPAN3.2; CD20L3; CDA01; MS4A6; MST090; MSTP090
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)