

Product datasheet for **GA112025**

Human MS4A5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MS4A5 CRISPRa kit - CRISPR gene activation of human membrane spanning 4-domains A5
Symbol:	MS4A5
UniProt ID:	<u>Q9H3V2</u>
Synonyms:	CD20-L2; CD20L2; TETM4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)