

Product datasheet for **GA112033**

Human ZNF323 (ZSCAN31) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZSCAN31 CRISPRa kit - CRISPR gene activation of human zinc finger and SCAN domain containing 31
Symbol:	ZSCAN31
UniProt ID:	Q96LW9
Synonyms:	dj874C20.2; FLJ23407; ZNF20-Lp; ZNF310P
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)