

Product datasheet for **GA112046**

Human HIF3 alpha (HIF3A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	HIF3A CRISPRa kit - CRISPR gene activation of human hypoxia inducible factor 3 subunit alpha
Symbol:	HIF3A
UniProt ID:	Q9Y2N7
Synonyms:	bHLHe17; HIF-3A; HIF3-alpha-1; IPAS; MOP7; PASD7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)