

Product datasheet for **GA112054**

Human **PAG608 (ZMAT3)** activation kit by **CRISPRa**

Product data:

| | |
|---------------|-------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits |
| Description: | ZMAT3 CRISPRa kit - CRISPR gene activation of human zinc finger matrin-type 3 |
| Symbol: | ZMAT3 |
| UniProt ID: | Q9HA38 |
| Synonyms: | PAG608; WIG-1; WIG1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)