


Product datasheet for **GA112057**

Human ZFP106 (ZNF106) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZNF106 CRISPRa kit - CRISPR gene activation of human zinc finger protein 106
Symbol:	ZNF106
Components:	GA112057G1, ZFP106 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112057G2, ZFP106 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112057G3, ZFP106 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H2Y7</u>
Synonyms:	SH3BP3; ZFP106; ZNF474
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)