

## Product datasheet for **GA112066**

### Human ARAP3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	ARAP3 CRISPRa kit - CRISPR gene activation of human ArfGAP with RhoGAP domain, ankyrin repeat and PH domain 3
Symbol:	ARAP3
UniProt ID:	<a href="#">Q8WWN8</a>
Synonyms:	CENTD3; DRAG1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)