


Product datasheet for **GA112076**

Human SDS3 (SUDS3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SUDS3 CRISPRa kit - CRISPR gene activation of human SDS3 homolog, SIN3A corepressor complex component
Symbol:	SUDS3
Components:	GA112076G1, SDS3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112076G2, SDS3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112076G3, SDS3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9H7L9
Synonyms:	SAP45; SDS3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)