

Product datasheet for **GA112077**

Human **TTC31** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TTC31 CRISPRa kit - CRISPR gene activation of human tetratricopeptide repeat domain 31
Symbol:	TTC31
Components:	GA112077G1, TTC31 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112077G2, TTC31 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112077G3, TTC31 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q49AM3</u>
Synonyms:	FLJ12788; FLJ33201; MGC120200
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)