

## Product datasheet for **GA112088**

### Human **TEKT3** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	TEKT3 CRISPRa kit - CRISPR gene activation of human tektin 3
Symbol:	TEKT3
UniProt ID:	<u><a href="#">Q9BXF9</a></u>
Synonyms:	FLJ32828
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)