

Product datasheet for **GA112106**

Human TMPRSS3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMPRSS3 CRISPRa kit - CRISPR gene activation of human transmembrane serine protease 3
Symbol:	TMPRSS3
UniProt ID:	P57727
Synonyms:	DFNB8; DFNB10; ECHOS1; TADG12
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)