

Product datasheet for **GA112107**

Human CSN7b (COPS7B) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	COPS7B CRISPRa kit - CRISPR gene activation of human COP9 signalosome subunit 7B
Symbol:	COPS7B
Components:	GA112107G1, CSN7b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112107G2, CSN7b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112107G3, CSN7b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9H9Q2
Synonyms:	CSN7B; SGN7b
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)