

Product datasheet for **GA112110**

Human PDI (PDIA2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PDIA2 CRISPRa kit - CRISPR gene activation of human protein disulfide isomerase family A member 2
Symbol:	PDIA2
UniProt ID:	Q13087
Synonyms:	PDA2; PDI; PDIP; PDIR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)