

Product datasheet for **GA112111**

Human RING finger protein unkempt like (UNKL) activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | UNKL CRISPRa kit - CRISPR gene activation of human unk like zinc finger |
| Symbol: | UNKL |
| UniProt ID: | Q9H9P5 |
| Synonyms: | C16orf28; ZC3H5L; ZC3HDC5L |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)