

Product datasheet for **GA112112**

Human WDR13 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR13 CRISPRa kit - CRISPR gene activation of human WD repeat domain 13
Symbol:	WDR13
Components:	GA112112G1, WDR13 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112112G2, WDR13 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112112G3, WDR13 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H1Z4</u>
Synonyms:	MG21
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)