

Product datasheet for **GA112117**

Human LPPR2 (PLPPR2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLPPR2 CRISPRa kit - CRISPR gene activation of human phospholipid phosphatase related 2
Symbol:	PLPPR2
UniProt ID:	Q96GM1
Synonyms:	LPPR2; PRG4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)