

Product datasheet for **GA112124**

Human Tensin 3 (TNS3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNS3 CRISPRa kit - CRISPR gene activation of human tensin 3
Symbol:	TNS3
Components:	GA112124G1, Tensin 3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112124G2, Tensin 3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112124G3, Tensin 3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q68CZ2</u>
Synonyms:	TEM6; TENS1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)