

Product datasheet for **GA112125**

Human RAI16 (FAM160B2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FAM160B2 CRISPRa kit - CRISPR gene activation of human family with sequence similarity 160 member B2
Symbol:	FHIP2B
Components:	GA112125G1, RAI16 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112125G2, RAI16 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112125G3, RAI16 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86V87</u>
Synonyms:	RAI16
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)