

## Product datasheet for **GA112130**

### Human S100P binding protein (S100BP) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	S100BP CRISPRa kit - CRISPR gene activation of human S100P binding protein
Symbol:	S100BPB
UniProt ID:	<u>Q96BU1</u>
Synonyms:	DKFZp313K2325; FLJ12903; S100PBPR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)