

## Product datasheet for **GA112160**

### Human P2Y12 (P2RY12) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	P2RY12 CRISPRa kit - CRISPR gene activation of human purinergic receptor P2Y12
Symbol:	P2RY12
UniProt ID:	<a href="#">Q9H244</a>
Synonyms:	ADPG-R; BDPLT8; HORK3; P2T(AC); P2Y(12)R; P2Y(AC); P2Y(ADP); P2Y(cyc); P2Y12; SP1999
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)