

Product datasheet for **GA112173**

Human **SLC13A3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	SLC13A3 CRISPRa kit – CRISPR gene activation of human solute carrier family 13 member 3
Symbol:	SLC13A3
Components:	GA112173G1, SLC13A3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112173G2, SLC13A3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112173G3, SLC13A3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8WWT9</u>
Synonyms:	NADC3; SDCT2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

