

Product datasheet for **GA112180**

Human PLEKHG2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLEKHG2 CRISPRa kit – CRISPR gene activation of human pleckstrin homology and RhoGEF domain containing G2
Symbol:	PLEKHG2
Components:	GA112180G1, PLEKHG2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112180G2, PLEKHG2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112180G3, PLEKHG2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H7P9</u>
Synonyms:	ARHGEF42; CLG; LDAMD
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

