

## Product datasheet for **GA112190**

### Human Lipin 3 (LPIN3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	LPIN3 CRISPRa kit - CRISPR gene activation of human lipin 3
Symbol:	LPIN3
UniProt ID:	<u><a href="#">Q9BQK8</a></u>
Synonyms:	LIPN3L; SMP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)