

Product datasheet for **GA112199**

Human **TTC23** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TTC23 CRISPRa kit - CRISPR gene activation of human tetratricopeptide repeat domain 23
Symbol:	TTC23
UniProt ID:	Q5W5X9
Synonyms:	1600012K10Rik; HCC-8, FLJ12572, FLJ12939; tetratricopeptide repeat domain 23
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)