

## Product datasheet for **GA112203**

### Human MRPS26 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | MRPS26 CRISPRa kit - CRISPR gene activation of human mitochondrial ribosomal protein S26                                                                                                                                                                                                                      |
| Symbol:       | MRPS26                                                                                                                                                                                                                                                                                                        |
| Components:   | GA112203G1, MRPS26 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA112203G2, MRPS26 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA112203G3, MRPS26 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q9BYN8</a></u>                                                                                                                                                                                                                                                                                 |
| Synonyms:     | C20orf193; dj534B8.3; GI008; MRP-S13; MRP-S26; MRPS13; NY-BR-87; RPMS13                                                                                                                                                                                                                                       |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)