

Product datasheet for **GA112205**

Human MRPS15 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MRPS15 CRISPRa kit - CRISPR gene activation of human mitochondrial ribosomal protein S15
Symbol:	MRPS15
Components:	GA112205G1, MRPS15 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112205G2, MRPS15 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112205G3, MRPS15 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P82914</u>
Synonyms:	DC37; MPR-S15; RPMS15; S15mt
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)