

Product datasheet for **GA112225**

Human P_{TOP} (ACD) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ACD CRISPRa kit - CRISPR gene activation of human ACD shelterin complex subunit and telomerase recruitment factor
Symbol:	ACD
UniProt ID:	Q96AP0
Synonyms:	PIP1; P _{TOP} ; TINT1; TPP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)