

Product datasheet for **GA112228**

Human ALS2CR4 (TMEM237) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM237 CRISPRa kit - CRISPR gene activation of human transmembrane protein 237
Symbol:	TMEM237
Components:	GA112228G1, ALS2CR4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112228G2, ALS2CR4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112228G3, ALS2CR4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q96Q45
Synonyms:	ALS2CR4; JBTS14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)