

Product datasheet for GA112231

Human MRPL44 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MRPL44 CRISPRa kit - CRISPR gene activation of human mitochondrial ribosomal protein L44
Symbol:	MRPL44
Components:	GA112231G1, MRPL44 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112231G2, MRPL44 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112231G3, MRPL44 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9H9J2
Synonyms:	COXPD16; L44MT; MRP-L44
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)