

Product datasheet for **GA112234**

Human TMEM135 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM135 CRISPRa kit - CRISPR gene activation of human transmembrane protein 135
Symbol:	TMEM135
Components:	GA112234G1, TMEM135 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112234G2, TMEM135 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112234G3, TMEM135 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86UB9</u>
Synonyms:	PMP52
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)