


Product datasheet for **GA112238**

Human UPF3B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	UPF3B CRISPRa kit - CRISPR gene activation of human UPF3B regulator of nonsense mediated mRNA decay
Symbol:	UPF3B
Components:	GA112238G1, UPF3B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112238G2, UPF3B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112238G3, UPF3B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BZ17</u>
Synonyms:	HUPF3B; MRX62; MRXS14; RENT3B; UPF3BP1; UPF3BP2; UPF3BP3; Upf3p-X; UPF3X
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)