

Product datasheet for **GA112239**

Human UPF3A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	UPF3A CRISPRa kit - CRISPR gene activation of human UPF3A regulator of nonsense mediated mRNA decay
Symbol:	UPF3A
Components:	GA112239G1, UPF3A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112239G2, UPF3A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112239G3, UPF3A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H1J1</u>
Synonyms:	HUPF3A; RENT3A; UPF3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)