

Product datasheet for **GA112243**

Human INTS3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	INTS3 CRISPRa kit - CRISPR gene activation of human integrator complex subunit 3
Symbol:	INTS3
Components:	GA112243G1, INTS3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112243G2, INTS3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112243G3, INTS3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q68E01</u>
Synonyms:	C1orf60; C1orf193; INT3; SOSS-A; SOSSA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)