

Product datasheet for **GA112244**

Human ANKRD57 (SOWAHC) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SOWAHC CRISPRa kit - CRISPR gene activation of human sosondowah ankyrin repeat domain family member C
Symbol:	SOWAHC
Components:	GA112244G1, ANKRD57 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112244G2, ANKRD57 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112244G3, ANKRD57 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q53LP3
Synonyms:	ANKRD57; C2orf26
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)