

Product datasheet for **GA112245**

Human WNK1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WNK1 CRISPRa kit - CRISPR gene activation of human WNK lysine deficient protein kinase 1
Symbol:	WNK1
UniProt ID:	Q9H4A3
Synonyms:	HSAN2; HSN2; KDP; p65; PPP1R167; PRKWNK1; PSK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)