

## Product datasheet for **GA112246**

### Human PCDH15 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | PCDH15 CRISPRa kit - CRISPR gene activation of human protocadherin related 15 |
| Symbol:       | PCDH15  |
| UniProt ID:   | <a href="#">Q96QU1</a>  |
| Synonyms:     | CDHR15; DFNB23; USH1F   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)       |

#### Product images:



[View online >](#)