

Product datasheet for **GA112250**

Human ZSWIM4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ZSWIM4 CRISPRa kit - CRISPR gene activation of human zinc finger SWIM-type containing 4
Symbol:	ZSWIM4
Components:	GA112250G1, ZSWIM4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112250G2, ZSWIM4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112250G3, ZSWIM4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H7M6</u>
Synonyms:	zinc finger, SWIM-type containing 4; zinc finger, SWIM domain containing 4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

