

Product datasheet for **GA112255**

Human PYCR3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PYCR3 CRISPRa kit - CRISPR gene activation of human pyrroline-5-carboxylate reductase 3
Symbol:	PYCR3
UniProt ID:	Q53H96
Synonyms:	PYCR3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)