


Product datasheet for **GA112257**

Human C8orf33 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C8orf33 CRISPRa kit - CRISPR gene activation of human chromosome 8 open reading frame 33
Symbol:	C8orf33
Components:	GA112257G1, C8orf33 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112257G2, C8orf33 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112257G3, C8orf33 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9H7E9
Synonyms:	FLJ20989
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)