

Product datasheet for **GA112262**

Human PLEKHA3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLEKHA3 CRISPRa kit - CRISPR gene activation of human pleckstrin homology domain containing A3
Symbol:	PLEKHA3
UniProt ID:	Q9HQB20
Synonyms:	FAPP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)