

Product datasheet for **GA112331**

Human TMEM38A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM38A CRISPRa kit - CRISPR gene activation of human transmembrane protein 38A
Symbol:	TMEM38A
Components:	GA112331G1, TMEM38A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112331G2, TMEM38A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112331G3, TMEM38A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9H6F2</u>
Synonyms:	TRIC-A; TRICA
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)