

Product datasheet for **GA112332**

Human TSEN34 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSEN34 CRISPRa kit - CRISPR gene activation of human tRNA splicing endonuclease subunit 34
Symbol:	TSEN34
Components:	GA112332G1, TSEN34 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA112332G2, TSEN34 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA112332G3, TSEN34 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9BSV6
Synonyms:	LENG5; PCH2C; SEN34; SEN34L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)