

Product datasheet for **GA112333**

Human KCTD15 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KCTD15 CRISPRa kit - CRISPR gene activation of human potassium channel tetramerization domain containing 15
Symbol:	KCTD15
UniProt ID:	Q96SI1
Synonyms:	MGC2628; MGC25497
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)