

Product datasheet for **GA112337**

Human TRPM8 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TRPM8 CRISPRa kit - CRISPR gene activation of human transient receptor potential cation channel subfamily M member 8
Symbol:	TRPM8
UniProt ID:	Q7Z2W7
Synonyms:	LTRPC6; TRPP8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)