

## Product datasheet for **GA112339**

### Human PRRG3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PRRG3 CRISPRa kit - CRISPR gene activation of human proline rich and Gla domain 3
Symbol:	PRRG3
UniProt ID:	<u><a href="#">Q9BZD7</a></u>
Synonyms:	PRGP3; TMG3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)